

#### **Features**

- Wide gain range: 1.5 to 200
- Independent integration and differentiation
- Complementary transconductance output
- ≤8 ns rise time (typically 6 ns)
- Ultra low noise
- Pole/Zero cancellation
- Preamp power output

### **Applications**

- Timing with Ge and Nal detectors
- General purpose wideband signal conditioning
- Signal to noise ratio optimization
- Ultra high count rate energy spectroscopy
- Rise time normalization

# **Model 2111** Timing Filter Amplifier

### Description

The Model 2111 Timing Filter Amplifier is a wideband, pulse shaping amplifier designed to enhance timing analysis. By adjusting the time constant settings, the user can optimize pulse shape and slope-to-noise ratio when making timing measurements. The Model 2111 can also serve as a wideband or variable bandwidth general purpose amplifier and as a shaping amplifier for high rate energy spectroscopy applications.

When used in a timing system, the Model 2111 will enhance performance of moderate (NaI) to slow (HPGe) speed radiation detectors. The  $\pm 5$  V output capability into a 50  $\Omega$ load makes the unit ideally suited to drive constant fraction discriminators such as the Model 2126. The selectable  $50/100 \Omega$  input impedance and a preamp power connector on the rear of the unit allow operation with a variety of preamps and detectors.

The gain of the Model 2111 spans the range from 1.5 to 200, using four coarse gain settings and a fine gain control for continuous coverage. The 2111's output stage is fully complementary and can drive both positive and negative unipolar pulses or bipolar pulses into 50  $\Omega$ . A front panel toggle switch is provided to select inverted or non-inverted output polarity. Independent integrate and differentiate controls allow flexible control of shaping with selection of time constants between 10 ns and 500 ns. With integrate and differentiate controls set in the Out position, wideband rise time is less than 8 ns with a Differentiate time constant greater than 150 µs. Pole/Zero cancellation overload recovery and high count rate performance for any Differentiate setting.

The flexible control of pulse shaping, the speed, and the wide dynamic range of the Model 2111 make it ideal for use in timing spectroscopy and other laboratory applications.

### **Specifications**

#### **INPUTS**

■ INPUT – Accepts positive or negative pulses or continuous signals from a detector, preamplifier, or other source; ac coupled;  $Z_{in}$  = 50 or 100  $\Omega$  for frequencies over 200 Hz, internal jumper selectable; dc input resistance 10 k $\Omega$ ; maximum input =  $\pm 5$  V divided by the amplifier gain, for linear output; protected to ±11 V ac, peak; 9 V ac, RMS; ±30 V dc.

#### **OUTPUTS**

■ OUTPUT – ±5 V into 50 Ω; short circuit proof; current output requires an external 50  $\Omega$  load.

## Model 2111 Timing Filter Amplifier

#### **CONTROLS**

- COARSE GAIN Rotary switch selects X3, X10, X30 or X100 gain.
- FINE GAIN Single-turn potentiometer; range X0.5 to X2.
- POLARITY Toggle switch selects INVerted or NON-INVerted OUTPUT polarity.
- DIFFERENTIATE Rotary switch selects OUT (no time constant), or 10, 20, 50, 100, 200 or 500 ns shaping; decay time in OUT position ≈ 160 µs.
- INTEGRATE Rotary switch selects OUT (no time constant), or 10, 20, 50, 100, 200 or 500 ns shaping; amplifier time constant in OUT position  $\approx 4$  ns.
- P/Z Twenty-turn screwdriver adjustable potentiometer cancels 20 µs to ∞ preamp time constants.

#### PERFORMANCE (with 50 $\Omega$ load)

- GAIN X1.5 to X200.
- RISE TIME ≤8 ns (typically 6 ns) with both shaping controls OUT or 2.2 times the selected INTEGRATE time constant.
- NOISE ≤10 μV with 200 ns INTEGRATE and DIFFERENTIATE time constants selected; ≤50 µV RMS, referred to input with filter OUT.
- BANDPASS 1 kHz to 45 MHz with both shaping controls OUT.
- NONLINEARITY ≤±1% integral.

- OVERSHOOT ≤10% with INTEGRATE OUT; ±2% with INTEGRATE time constant selected.
- GAIN DRIFT <0.1%/°C.
- DC DRIFT <25 µV/°C.</p>
- OUTPUT OFFSET ≤3 mV.

#### **CONNECTORS**

- SIGNAL INPUT and OUTPUT connectors are BNC type.
- PREAMP POWER Amphenol type 17-10070; rear panel.

#### **POWER REQUIREMENTS**

+24 V dc - 55 mA +12 V dc - 170 mA\* -24 V dc - 80 mA

-12 V dc - 150 mA

\*This exceeds the normal Bin power allotment of 167 mA for a single width module.

- SIZE Standard single-width NIM module; 3.43 x 22.12 cm (1.35 x 8.71 in.) per DOE/ER-0457T.
- NET WEIGHT 0.9 kg (2.0 lb).
- SHIPPING WEIGHT 1.8 kg (4.0 lb).

#### **ENVIRONMENTAL**

- OPERATING TEMPERATURE 0 to 50 °C.
- OPERATING HUMIDITY 0 to 80% relative, non-condensing. Meets the environmental conditions specified by EN 61010, Installation Category I, Pollution Degree 2.




Figure 1 Timing with an HPGe Detector

# Model 2111 Timing Filter Amplifier




Figure 2a Differentiation with Integration OUT Differentiation set at OUT, 10, 20, 50, 100, 200 and 500 ns




Figure 2c **Identical Integration and Differentiation settings** Controls set at OUT, 10, 20, 50, 100, 200 and 500 ns




Figure 2b Integration with Differentiation OUT Integration set at OUT, 10, 20, 50, 100, 200 and 500 ns





